

Self-developed control intelligent system, it can manage various functions in real time. The equipment allows for precise integration, providing a more convenient operating experience and enhancing efficiency.

Stability

Timotion's imported motor ensures seamless lifting, operating silently and stably without jamming, which provides patients with a sense of safety during medical treatment.

Weighty

Integrated precision steel bending plate, with strong bearing capacity, heavy weight, and strong support, can withstand even the largest weight

Convenient

After extensive research on ergonomics, we have designed the most suitable dental chair for both dentist and patients.

Durable

The dental chair is constructed from high-quality materials sourced both domestically and internationally, ensuring durability and relevance over time, with no concerns about after-sales service.



Zinc silicone strong and weak suction

Allows 135 C sterilization, Aging-resistant silicone, impact resistance



Stainless steel interface

Dust-proof and tear-resistant tube



Galvanized valve

High precision, rust and corrosion resistance, long service life, low maintenance rate



Rotary Assistant Tray

Four hands operation, easy control



Push-button armrests

Prevent accidental touch, safer dentist-patient interaction



Deluxe Ceramic detachable spittoon

90°rotary, easy to clean, reduce risk of cross-infection

PU Color options »



H01

Light grey



H02

Dark grey



H03

Middle blue



H04

Dark blue



H05

Apple green



H06

Purplish red



H07

Rose red



H08

Yellow



H09

Orange



H10

Bright Yellow

**Skilled in stability
Make quality visible**

The chair frame design features a precise structural layout. After conducting numerous tests, we have refined the manufacturing process, establishing a stable support system that feels as solid and enduring as Mount Tai. The operation and controls are smoother, providing patients with enhanced comfort both physically and emotionally.



C300

Iron hand in the velvet glove

Experience beauty inside and out



Comfortable dental unit

C300

Smart Offers for Saving and Sharing

New Upgrade | Increased Performance
Complete Configuration | Highly Competitive



Digital Sensor Oral lamp

Philips LED chip | Injection Plastic Cover
Backcast Aluminum Plate



Two Adjustment Methods:
Illuminance and Color temperature

Illuminance:
Adjustable from 6000Lux-40000Lux

Color Temperature:
Adjustable from 5000k-5700k


Brightness Memory Function:
When the Oral Lamp is turned off after normal use, it will automatically adjust the brightness from the lowest setting to the brightness mode set before the last shutdown the next time it is turned on.

Adjustable height range of dental chair

Stable Running | Smoothly Moving

Gentle Movement

Start/Stop
Integrated with Rostar's self-developed soft start system, it eliminates frustration when the motor starts and stops, allowing the dental chair to flow smoothly and freely. This system also reduces the patient's discomfort and tension during treatment.



Backrest running smoothly and efficiently

Ingenuity design



The backrest has excellent resetting performance and is easy to use. With one button, dentists can quickly adjust it to the lowest or highest position they need. This helps them work more efficiently and makes the patient more comfortable during the entire process.

Considerate and Caring

Ergonomic Design | Strong Support
All-round Treatment



Super Resistant PU & High Elastic Seat Cushion

Technology Fabric | Anti-mold And Antibacterial
High Elastic Inner Core | Not Easy To Crack

Extreme Ergonomics Described in detail

The extra-large instrument tray panel design allows dentists to operate more flexibly and conveniently. Silicone pads help prevent dirt and bacteria accumulation.

- Extra large storage capacity
- Ultra-Stable Panel
- Anti-slip and dirt-resistant silicone pad
- Smooth anti-collision corners



Extended seat cushion More compatible

Ergonomic design | Human Height coverage
Extended seat cushions accommodate various heights, providing more relaxed and comfortable support for the legs.



Super resistant PU leather clean and hygienic

Firm fit | Easy to clean

Simply wipe it with a rag dampened with water to reduce the time needed for medical care and maintenance, while extending the service life of the chair.