

# C200



## Smart Offers for Saving and Sharing

New Upgrade | Increased Performance  
Complete Configuration | Highly Competitive

### LED Sensor Oral lamp

Large Spot | Strong brightness  
Soft Light | Protect Eyes

Two Adjustment Methods  
Infinitely Adjustable Brightness  
Full Coverage Lighting  
Maximum Rotation 90°



### Adjustable height range of dental chair

Stable Running | Smoothly Moving

#### Gentle Movement

##### Start/Stop

Integrated with Rostar's self-developed soft start system, it eliminates frustration when the motor starts and stops, allowing the dental chair to flow smoothly and freely. This system also reduces the patient's discomfort and tension during treatment.



Lowest chair position  
> 390mm

Highest chair position  
< 680mm

### Backrest running smoothly and efficiently

Ingenuity design



The backrest has excellent resetting performance and is easy to use. With one button, dentists can quickly adjust it to the lowest or highest position they need. This helps them work more efficiently and makes the patient more comfortable during the entire process.



### Considerate and Caring

Ergonomic Design | Strong Support  
All-round Treatment

### Super Resistant PU & High Elastic Seat Cushion

Technology Fabric | Anti-mold And Antibacterial  
High Elastic Inner Core | Not Easy To Crack

### Extreme Ergonomics Described in detail

The extra-large instrument tray panel design allows dentists to operate more flexibly and conveniently. Silicone pads help prevent dirt and bacteria accumulation.

- Extra large storage capacity
- Ultra-Stable Panel
- Anti-slip and dirt-resistant silicone pad
- Smooth anti-collision corners



### Extended seat cushion More compatible

Ergonomic design | Human Height coverage  
Extended seat cushions accommodate various heights, providing more relaxed and comfortable support for the legs.

### Super resistant PU leather clean and hygienid

Firm fit | Easy to clean  
Simply wipe it with a rag dampened with water to reduce the time needed for medical care and maintenance, while extending the service life of the chair.



### Zinc silicone strong and weak suction

Allows 135°C sterilization, Aging-resistant silicone, impact resistance



### Rotation unit box

Increase the treatment space for convenient four-handed operation




### New Upgrade Multifunctional pedal

New upgrade with complete functionality.



### Stainless steel interface

Dust-proof and tear-resistant tube




### Galvanized valve

High precision, rust and corrosion resistance, long service life, low maintenance rate



### Instrument Tray

Clear and concise keys that are elegant in design. Injection durable plastic




### Rotary Assistant Tray

Four hands operation, easy control



### Rotatable handrail

It is convenient and considerate to meet patients' needs for getting on and off the chair



### Instrument Tray Pipe diagram

Water and Air pipe imported from BASF, USA, with lifetime warranty

### PU Color options »



# C200

Activate Your Quality Experience  
Empowering Your Treatment

Comfortable dental unit

